

| | |
|--------------------------------------|-----------|
| Guide Wires | 6.1 - 6.2 |
| Single-use E.R.C.P. Catheters | 6.3 |
| Reusable E.R.C.P. Catheters | 6.4 - 6.6 |

Guide Wires &
E.R.C.P. Catheters





STERILE EO

Nitinol Guide Wire "Hi-Flex"

- straight tip
- 5 cm hydrophilic tip



| Order No. | diameter | length | flexible part | colour |
|-----------|-----------|--------|---------------|----------------|
| 217 35450 | .035 inch | 450 cm | 5 cm | violet / white |

Nitinol Guide Wire "Hi-Flex"

- curved tip
- 5 cm hydrophilic tip



| Order No. | diameter | length | flexible part | colour |
|-----------|-----------|--------|---------------|----------------|
| 218 35450 | .035 inch | 450 cm | 5 cm | violet / white |

Nitinol Guide Wire "Smart-Guide"

- 50 cm hydrophilic coating



| Order No. | diameter | length | flexible part | colour |
|-----------|-----------|--------|---------------|----------------|
| 216 35450 | .035 inch | 450 cm | 5 cm | violet / white |

134°C
273 F


STERILE EO

Nitinol Guide Wires "Stripe-Guide"

- reusable & autoclavable



| Order No. | diameter | length | flexible part | colour |
|-----------|-----------|--------|---------------|-----------------|
| 215 25400 | .025 inch | 400 cm | 5 cm | violet / yellow |
| 215 35400 | .035 inch | 400 cm | 5 cm | violet / white |
| 215 35500 | .035 inch | 500 cm | 5 cm | violet / white |





STERILE EO

Guide Wire TERUMO® RADIOFOCUS®

- straight tip



| Order No. | diameter | length | flexible part |
|---------------|-----------|--------|---------------|
| 900 25 263GSM | .025 inch | 260 cm | 3 cm |
| 900 35 263GSM | .035 inch | 260 cm | 3 cm |
| 900 35 403GSM | .035 inch | 400 cm | 3 cm |

Guide Wire TERUMO® RADIOFOCUS®

- curved tip



| Order No. | diameter | length | flexible part |
|---------------|-----------|--------|---------------|
| 900 25 263GAM | .025 inch | 260 cm | 3 cm |
| 900 35 263GAM | .035 inch | 260 cm | 3 cm |
| 900 35 403GAM | .035 inch | 400 cm | 3 cm |


Ø 1.8
mm



STERILE EO

E.R.C.P. Catheter, standard tip

- accepts guide wire .035 inch
- radiopaque tip



packing unit

5 pieces

E 43 18-GW

E.R.C.P. Catheter, tapered tip

- accepts guide wire .035 inch
- radiopaque tip



packing unit

5 pieces

E 44 18-GW

E.R.C.P. Catheter, hour glass metal tip

- accepts guide wire .035 inch
- radiopaque tip



packing unit

5 pieces

E 48 18-GW



single-use


E.R.C.P. Catheters

Ø 1.8 mm

134 °C
273 °F

E.R.C.P. Catheter, standard tip

- accepts guide wire .035 inch



with mandrin

with mandrin / proximal Ø 2.45 mm

43 18

43 22-18FD



E.R.C.P. Catheter, tapered tip

- accepts guide wire .035 inch



with mandrin

with mandrin / proximal Ø 2.45 mm

44 18

44 22-18FD



E.R.C.P. Catheter, conical metal tip

- accepts guide wire .035 inch



with mandrin

with mandrin / proximal Ø 2.45 mm

45 18

45 22-18FD



Ø 1.8 mm

134 °C
273 °F

E.R.C.P. Catheter, tapered metal tip

- accepts guide wire .035 inch



with mandrin

with mandrin / proximal Ø 2.45 mm

46 18

46 22-18FD



E.R.C.P. Catheter, round metal tip

- accepts guide wire .035 inch



with mandrin

with mandrin / proximal Ø 2.45 mm

47 18

47 22-18FD



E.R.C.P. Catheter, hour-glass metal tip

- accepts guide wire .035 inch



with mandrin

with mandrin / proximal Ø 2.45 mm

48 18

48 22-18FD



Ø 1.8 mm

134°C
273 F

needle tip 2 mm

- not to be used with guide wire



with mandrin

49 18



Ø 1.6 mm

134°C
273 F

tapered tip

- accepts guide wire
.025 inch



with mandrin

44 16

