

Single-Use	<b>Esophageal Multiband Ligator “MBL”</b>	12.1
Reusable	<b>Suction Ligation Device</b>	12.2
Single-Use	<b>Clip-Applicators</b>	12.3
Reusable	<b>Spray Catheters</b>	12.4

Multiband Ligators, Clip-  
Applicators, Spray Catheters





STERILE EO

### Esophageal Multiband Ligator "MBL"

- pre-loaded barrel with 6 rings
- improved mounting to the endoscope
- additional luer-lock irrigation port
- suitable for all endoscopes from 9.5 to 11.5 outer diameter

Order No.

MBL-01-06



single-use

Esophageal Multiband Ligator "MBL"

134 °C  
273 °F

### Suction Ligation Device for ligation of haemorrhoids 2. and 3. degree

- including:
- large suction head
- small suction head
- 2 loading cones
- 100 ligation rings

Order No.

178 000	Suction Ligation Device, complete set
---------	---------------------------------------

### Individual Components

Order No.

178 102	loading cone
178 110	small suction head
178 111	large suction head
178 200	ligation rings (1 package = 100 pieces)






STERILE EO

### Clip-Applicator

- Single-use Clip-Applicator for endoscopic marking & hemostasis
- rotatable & repositionable
- re-opening capability

Order No.	length	jaw span	arm length	for working channel
CLP-12-135-B	160 cm	12.5 mm	5 mm	2.8 mm and larger
CLP-12-135-G	230 cm	12.5 mm	5 mm	2.8 mm and larger



single-use

Clip-Applicators

Ø 2.3  
mm

134 °C  
273 °F

#### radial Spray Catheter

- for marking suspicious mucosa
- 360° radial spraying
- LL connector with mandrin

length	Order No.
120 cm	47 226-A
160 cm	47 226-B
180 cm	47 226-C
230 cm	47 226-G



Ø 1.8  
mm

134 °C  
273 °F

#### Spray Catheter for wide spray

- for marking suspicious mucosa
- for straight, cone-shaped spraying
- LL connector with mandrin

length	Order No.
120 cm	47 185-A
160 cm	47 185-B
180 cm	47 185-C
230 cm	47 185-G

