

# 60+, 100+ & 200+ fingerswitches

## 3m

**PMS 550**  
200 Plus fingerswitch  
with blade electrode  
pack 1

plus  
**AUTOClabel:**

**PMS 750**  
100 Plus fingerswitch  
with blade electrode  
pack 5

plus  
**AUTOClabel:**

**PMS 950**  
60 Plus fingerswitch  
with blade electrode  
pack 5

plus  
**AUTOClabel:**

## 5m

**PMS 550/5**  
200 Plus fingerswitch  
with blade electrode  
pack 1

plus  
**AUTOClabel:**

**PMS 750/5**  
100 Plus fingerswitch  
with blade electrode  
pack 5

plus  
**AUTOClabel:**

**PMS 950/5**  
60 Plus fingerswitch  
with blade electrode  
pack 5

plus  
**AUTOClabel:**

## 200+



Autoclavable for a maximum of 200 cycles

## 100+



Autoclavable for a maximum of 100 cycles

## 60+



Autoclavable for a maximum of 60 cycles

## AUTOClabels




Autoclave tracking system



cable terminated  
with 3-pin plug

for adaptors see page 7




Ergonomically designed for precise control



Non slip grips - for safety



**AUTOCLABELS** autoclave tracking system



Non slip grips - for safety

Produced in environmentally friendly polypropylene (PP)

Universal connector fits all 2.38mm and 4.00mm active electrodes (USA and European fittings)

### Ordering Details:

Pack of 3 footswitch controlled handles with blade electrode, AUTOCLABEL and punch, terminated with electrosurgery generator connector:

**PMS 651**  
with 3m cable



Terminated with connector to fit:  
**4mm** Aesculap  
Berchtold  
Erbe T. Series  
Eschmann TD300, TD311, DS302  
Sutter BM780

**PMS 651/5**  
with 5m cable

**PMS 652**  
with 3m cable



Terminated with connector to fit:  
**8mm** Bard, Birtcher  
Bovie, Bowa  
Conmed, Eschmann  
Medtrex  
Valleylab

**PMS 652/5**  
with 5m cable

**PMS 654**  
with 3m cable



Terminated with connector to fit:  
**5mm** ACC (standard)  
ICC  
Erbe V10\*  
Karl Storz  
\*Depends on generator socket

**PMS 654/5**  
with 5m cable

**PMS 653**  
with 3m cable



Terminated with connector to fit:  
**4mm (female)**  
Aesculap,  
Berchtold  
Erbe T Series,  
Huttlinger  
Martin, Micromed  
Siemens  
Sutter BM780

**PMS 653/5**  
with 5m cable

# Fingerswitch adaptors

## 3-pin plug to Berchtold / Martin generators

### PMS 800

Pack of 1  
Fingerswitch adaptor



## 3-pin plug to Erbe ACC/ICC generators with 5mm pin

### PMS 850

Pack of 1  
Fingerswitch adaptor



## 3-pin plug to Erbe T series generators with 4mm pin

### PMS 860

Pack of 1  
Fingerswitch adaptor



## Connector from fingerswitch



cable terminated  
with 3-pin plug